

Recombinant Human Uridine Phosphorylase 1,UPP1 (N-6His)

Catalog No: #AP75814

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75814-1 10ug #AP75814-2 50ug #AP75814-3 500ug #AP75814-4 1mg

Description

Product Name	Recombinant Human Uridine Phosphorylase 1,UPP1 (N-6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q86Y75
Uniprot	Q86Y75
Calculated MW	21.3kDa
Target Sequence	MGSSHHHHHSSGLVPRGSHMQRKLVTSLEPGTVVITEQAVDTCFKAEFEQIVLGKRVIRKTDLNKKLVQEL LLCSAELSEFTTVVGNTMCTLDFYEGQGRLDGLCSYTEKDKQAYLEAAYAAGVRNIEMESSVFAAMCSACG LQAAVVCVTLNRLLEGDQISSPRNVLSEYQQRQPQLVSYFIKKLSKA
Formulation	Supplied as a 0.2 µm filtered solution of 20mM Tris, 200mM NaCl, 1mM DTT, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only