

Recombinant Human C1QBP,HABP1 (C-6His)

Catalog No: #AP75815



Package Size: #AP75815-1 10ug #AP75815-2 50ug #AP75815-3 500ug #AP75815-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human C1QBP,HABP1 (C-6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q07021
Uniprot	Q07021
GeneID	708;
Calculated MW	24.9kDa
Target Sequence	MLHTDGDKAFVDFLSDEIKEERKIQKHKTLPKMSSGGWELELNGTEAKLVRKVAGEKITVTFNINNSIPPTFDGE EEPSQGQKVVEEQPELTSTPNFVVEVIKNDGKKALVLDCHYPEDEVGQEDEAESDIFSIREVVSFQSTGESEW KDTNYTLNTDSLWALYDHLMDFLADRGVDNTFADELVELSTALEHQEYITFLEDLKSFKVKSQLEHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM Tris, 20% Glycerol, 1mM DTT, pH 7.5.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only