

Recombinant Human Sulfotransferase 4A1,SULT4A1

Catalog No: #AP75816



Package Size: #AP75816-1 10ug #AP75816-2 50ug #AP75816-3 500ug #AP75816-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Sulfotransferase 4A1,SULT4A1
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9BR01
Uniprot	Q9BR01
GeneID	25830;
Calculated MW	33kDa
Target Sequence	MAESEAETPSTPGEFESKYFEFHGVRLPPFCRGKMEEIANFPVRPSDVWIVTYPKSGTSSLQEVVYLVSQGAD PDEIGLMNIDEQLPVLEYPQPGLDIIKELTSPRLIKSHLPYRFLPSDLHNGDSKVIYMARNPKDLVVSYYQFHRS RTMSYRGTFQEFRRFMNDKLGYSWFEHVQEFWEHRMDSNVFLKYEDMHRDLVTMVEQLARFLGVSCD KAQLEALTEHCHQLVDQCCSAEALPVGRGRVGLWKDIFTVSMNEKFDLVYKQKMGKCDLTFDFYL
Formulation	Supplied as a 0.2 µm filtered solution of 20mM PB,150mM NaCl,pH7.4.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only