

Recombinant Human DnaJ Homolog Subfamily B Member 1, DNAJB1 (C-6His)



Catalog No: #AP75817

Orders: order@signalwayantibody.com

Package Size: #AP75817-1 10ug #AP75817-2 50ug #AP75817-3 500ug #AP75817-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human DnaJ Homolog Subfamily B Member 1, DNAJB1 (C-6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P25685
Uniprot	P25685
GeneID	3337;
Calculated MW	39.1kDa
Target Sequence	GKDYYQTLGLARGASDEEIKRAYRRQALRYHPDKNKEPGAEEKFKEIAEAYDVLSDPRKREIFDRYGEEGLKG SGPSGGSGGGANGTSFSYTFHGDPHAMFAEFFGGRNPFDTFFGQRNGEEMDIDDPFSGFPMGMGGFTNV NFGRSRSQAQEPARKKQDPPVTHDLRVSLLEIYSGCTKKMKISHKRLNPDGKSIRNEDKILTIEVKKGWKEGTKI TFPKEGDQTSNNIPADIVFLKDKPHNIFKRDGSDVIYPARISLREALCGCTVNVPTLDGRTIPVVKDVIKPGM RRKVPGEGLPLPKTPEKRGDLIEFEVIFPERIPQTSRTVLEQVLPILHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, 1mM EDTA, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only