

# Recombinant Human Synaptobrevin Homolog YKT6, YKT6 (N-6His)

Catalog No: #AP75818

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP75818-1 10ug #AP75818-2 50ug #AP75818-3 500ug #AP75818-4 1mg

## Description

Product Name	Recombinant Human Synaptobrevin Homolog YKT6, YKT6 (N-6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	O15498
Uniprot	O15498
GeneID	10652;
Calculated MW	24.5kDa
Target Sequence	MNHKVVHHHHHMKLYSLSVLYKGEAKVLLKAAAYDVSSFSFFQRSSVQEFMTFTSQLIVERSSKGTRASVKE QDYLCHVYVRNDSLAVVIADNEYPSRVAFTLLEKVLDEFKQVDRIDWPVGSPTIHYPALDGHLSTRYQNPR EADPMTKVQAELDETKIILHNTMESLLERGEKLDLVLKSEVLGTQSKAFYKTARKQNSCCAIMLDLQSR
Formulation	Supplied as a 0.2 µm filtered solution of 50mM Tris, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only