

Recombinant Human Actin-Related Protein 8,ACTR8 (N-6His)



Catalog No: #AP75821

Orders: order@signalwayantibody.com

Package Size: #AP75821-1 10ug #AP75821-2 50ug #AP75821-3 500ug #AP75821-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Actin-Related Protein 8,ACTR8 (N-6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9H981-3
Uniprot	Q9H981
GeneID	93973;
Calculated MW	38.1kDa
Target Sequence	MNHKVVHHHHMILMKMGFSGIVVHQESVCATYGSGLSSTCIVDVGDKQTSVCCVEDGVSHRNTRIFSWNQ DISGLQDHEFQIRHPDSPALLYQFRLGDEKLQAPMALFYPATFGIVGQKMTTLQHRSQGDPEDPHDEHYLLAT QSKQEQSAKATADRKSASKPIGFEGDLRGQSSDLPERLHSQEVDLGSAAQGDGLMAGNDSEEALTALMSRKT AISLFEGKALGLDKAILHSIDCCSSDDTKKKMYSSILVVGGGLMFHKAQEFLQHRILNKMPFSFRRIENVVITR PKDMDPRLIAWKGGAVLAFLDQELWIYQREWQRFGVRMLRERAAFVW
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB,150mM NaCl,1mM DTT,pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only