

Recombinant Human Sulfotransferase 1B1,SULT1B1 (N-6His)

Catalog No: #AP75822

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75822-1 10ug #AP75822-2 50ug #AP75822-3 500ug #AP75822-4 1mg

Description

Product Name	Recombinant Human Sulfotransferase 1B1,SULT1B1 (N-6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	O43704
Uniprot	O43704
GeneID	27284;
Calculated MW	36.3kDa
Target Sequence	MNHKVVHHHHHMLSPKDILRKDLKLVHGYPMTCAFASNWEKIEQFHSRPDDIVIATYPKSGTTWVSEIIDMILN DGDIEKCKRGFITEKVPMLEMTLPGLRTSGIEQLEKNPSRIVKTHLPTDLLPKSFWENNCKMIYLARNAKDVS VSYHFDFLMNNLQFPFGTWEEYLEKFLTGKVAYGSWFTHVKNWWKRKEEHPIFLYEDMKENPKKEIKKIIR FLEKNLNDEILDRIIHHTSFEVMKDNPLVNYTHLPTTVMHDSKSPFMRKGTAGDWKNYFTVAQNEKFDIAYETE MSKTALQFRTEI
Formulation	Supplied as a 0.2 µm filtered solution of 20mM Tris, 0.1M NaCl, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only