

# Recombinant Human Sulfotransferase 2A1,SULT2A1 (N-6His)

Catalog No: #AP75824

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP75824-1 10ug #AP75824-2 50ug #AP75824-3 500ug #AP75824-4 1mg

## Description

Product Name	Recombinant Human Sulfotransferase 2A1,SULT2A1 (N-6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q06520
Uniprot	Q06520
GeneID	6822;
Calculated MW	35.2kDa
Target Sequence	MNHKVVHHHHHMSDDFLWFEGIAFPTMGFRSETLRKVRDEFVIRDEDVILTYPKSGTNWLAEILCLMHSKGD AKWIQSVPIWERSPWVESEIGYTALSETESPRLFSSHLPIQLFPKSFSSKAKVIYLMRNPRDVLVSGYFFWKN MKFIKKPKSWEEYFEWFCQGTVLYGSWFDHIHGWMMPMREEKNLLLLSYEELKQDTGRTIEKICQFLGKLEPE ELNLILKNSSFQSMKENKMSNYLLSVDYVVDKAQLLRKGVSGDWKNHFTVAQAEDFDKLFQEKMADLPREL FPWE
Formulation	Supplied as a 0.2 µm filtered solution of 20mM Tris, 100mM NaCl, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only