

## Recombinant Human Caspase-10,CASP10 (C-6His)

Catalog No: #AP75825



Package Size: #AP75825-1 10ug #AP75825-2 50ug #AP75825-3 500ug #AP75825-4 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                 |  |
|-----------------|--|
| Product Name    | Recombinant Human Caspase-10,CASP10 (C-6His)   |
| Host Species    | E. coli  |
| Purification    | Greater than 95% as determined by reducing SDS-PAGE.   |
| Specificity     | Human  |
| Accession No.   | Q92851-2   |
| Uniprot         | Q92851   |
| GeneID          | 843;   |
| Calculated MW   | 30.1kDa  |
| Target Sequence | MVKTFLEALPQESWQNKHAGSNGNRATNGAPSLVSRGMQGASANTLNSETSTKRAAVYRMNRNHRGLCVIV<br>NNHSFTSLKDRQGTHKDAEILSHVFQWLGFTVHIHNNVTKVEMEMVLQKQKCNPAHADGDCFVFCILTHGRF<br>GAVYSSDEALIPREIMSHFTALQCPRLAEKPKLFFIQACQGEEIQPSVSIADALNPEQAPTSLQDSIPAEADFL<br>LGLATVPGYVSFRHVVEGSWYIQSLCNHLKLVPRHEDILSILEHHHHHH |
| Formulation     | Supplied as a 0.2 µm filtered solution of 25mM HEPES, 10mM DTT, pH 7.5.  |
| Storage         | Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.   |

Note: This product is for in vitro research use only