

Recombinant Human Serpin B5, SERPINB5, Maspin

Catalog No: #AP75826



Package Size: #AP75826-1 10ug #AP75826-2 50ug #AP75826-3 500ug #AP75826-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Serpin B5, SERPINB5, Maspin
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P36952-2
Uniprot	P36952
GeneID	5268;
Calculated MW	25.7kDa
Target Sequence	MDALQLANSFAVDLQKQCEKEPLGNVLFSPICLSTSLSLAQVGAKGDTANEIGQVLHFENVKDVPFQFQTVT SDVNKLSFYSLKLIKRLYVDKSLNLSTEFISSTKRPYAKELETVDKDKLEETKQGQINNSIKDLTDGHFENILAD NSVNDQTKILVVNAAYFVGKWMKKFPESETKECPFRVNVKCGAACSSKRSPIIDVKNDRDRVGHKSIPMRNLR ARPAKCLS
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only