

Recombinant Human Sulfotransferase 1A1,SULT1A1 (N-6His)

Catalog No: #AP75828

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75828-1 10ug #AP75828-2 50ug #AP75828-3 500ug #AP75828-4 1mg

Description

Product Name	Recombinant Human Sulfotransferase 1A1,SULT1A1 (N-6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P50225
Uniprot	P50225
GeneID	6817;
Calculated MW	35.5kDa
Target Sequence	MNHKVVHHHHHMELIQDTSRPPLEYVKGVPLIKYFAEALGPLQSFQARPDLLISTYPKSGTTWVSQILDMIYQ GGDLEKCHRAPIFMRVPFLEFKAPGIPSGMETLKDTPAPRLLKTHLPLALLPQTLLDQKVVVYVARNADVAV SYYHFYHMAKVHPEPGTWDSFLEKFMVGEVSYGSWYQHVQEWWELSRTHPVLYLFYEDMKENPKREIQKIL EFVGHSLPEETVDFVVQHTSFKEMKKNPMTNYTTVPQEFMDHSISPFRKGMAGDWKTTFTVAQNERFDAD YAEKMAGCSLSFRSEL
Formulation	Supplied as a 0.2 µm filtered solution of 20mM Tris,150mM NaCl,1mM EDTA,pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only