

Recombinant Human N-Acetylneuraminase Synthase, NANS, SAS (N-6His)

Catalog No: #AP75829

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75829-1 10ug #AP75829-2 50ug #AP75829-3 500ug #AP75829-4 1mg

Description

Product Name	Recombinant Human N-Acetylneuraminase Synthase, NANS, SAS (N-6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9NR45
Uniprot	Q9NR45
GeneID	54187;
Calculated MW	42.4kDa
Target Sequence	MGSSHHHHHHSSGLVPRGSHMPLELELCPGRWVGGQHPCFIIAEIGQNHQGDLDAKRMIRMAKECGADCA KFQKSELEFKFNRKALDRPYTSKHSWGKTYGEHKRHLEFSDQYRELQRYAEEVGIFFTASGMDEMAVEFLH ELNVPFFKVGSGDTNFPYLEKTAKKGRPMVISSGMQSMQVYQIVKPLNPNFCFLQCTSAYPLQPEDV NLRVISEYQKLFDPDIPIGYSGHETGIAISVAVALGAKVLERHITLTKTWKGDHSASLEPGELAEVRSVRLVER ALGSPTKQLLPCEMACNEKLGKSVVAKVKIPEGTILTMMLTVKVGEPKGYPPEDIFNLVGKKVLVTVVEEDDTI MEELVDNHGKKIKS
Formulation	Supplied as a 0.2 µm filtered solution of 20mM Tris, 100mM NaCl, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only