

# Recombinant Human Ubiquitin-Conjugating Enzyme E2 C,UBE2C,UBCH10 (N-6His)



Catalog No: #AP75831

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP75831-1 10ug #AP75831-2 50ug #AP75831-3 500ug #AP75831-4 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human Ubiquitin-Conjugating Enzyme E2 C,UBE2C,UBCH10 (N-6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	O00762
Uniprot	O00762
GeneID	11065;
Calculated MW	23.3kDa
Target Sequence	MRGSHHHHHHGSMASQNRDPAATSVAAARKGAEPSSGAARGPVGKRLQQELMTLMMMSGDKGISAFPESDN LFKVVGTIHGAAGTVYEDLRYKLSLEFSPGYPNAPTVMKFLTPCYHPNVDTQGNICLDILKEKWSALYDVRTILL SIQSLLGEPNIDSPLNTHAAELWKNPTAFKKYLQETYSKQVTSQEP
Formulation	Supplied as a 0.2 µm filtered solution of 20mM PB,150mM NaCl, pH 7.0.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only