

Recombinant Human Amphiregulin,AREG

Catalog No: #AP75833



Package Size: #AP75833-1 10ug #AP75833-2 50ug #AP75833-3 500ug #AP75833-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Amphiregulin,AREG
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P15514
Uniprot	P15514
GeneID	374;
Calculated MW	11.4kDa
Target Sequence	MSVRVEQVVKPPQNKTESENTSDKPKRKKKGGKNGKNRRNRKKKNPCNAEFQNFICIHGECKYIEHLEAVTC KCQQEYFGERCGEKSMKTHSMIDSSLSK
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB,150mM NaCl,pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only