

Recombinant Human PDCD5,TFAR19 (N-6His)

Catalog No: #AP75834



Package Size: #AP75834-1 10ug #AP75834-2 50ug #AP75834-3 500ug #AP75834-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human PDCD5,TFAR19 (N-6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	O14737
Uniprot	O14737
GeneID	9141;
Calculated MW	16.4kDa
Target Sequence	MGSSHHHHHHSSGLVPRGSHMADEELEALRRQRLAELQAKHGDPGDAAQQAQKHEAEMRNSILAQVLDQ SARARLSNLALVKPEKTKAVENYLIQMARYGQLSEKVSEQGLIEILKVSQQTEKTTTVKFNRRKVMDSDEDD DY
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB150mM NaCl,pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only