

Recombinant Human Myozenin-2, MYOZ2 (C-6His)

Catalog No: #AP75835



Package Size: #AP75835-1 10ug #AP75835-2 50ug #AP75835-3 500ug #AP75835-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Myozenin-2, MYOZ2 (C-6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9NPC6
Uniprot	Q9NPC6
GeneID	51778;
Calculated MW	30.9kDa
Target Sequence	<p> MLSHNTMMKQRKQQATAIMKEVHGNDVDGMDLGKKVSIPRDIMLEELSHLSNRGARLFKMRQRRSDKYTFE NFQYQSRAQINHSIAMQNGKVDGSNLEGGSQQAPLTPPNTDPDRSPPNPDNIAPGYSGPLKEIPPEKFNTTAV PKYYQSPWEQAIENDPELLEALYPKLFKPEGKAELPDYRSFNRVATPFGGFEEKASRMVFKVPDFELLLLDTP RFMSFVNPLSGRRSFNRTPKGWISENIPVITTEPTDDTTVPESEDLLEHHHHHH </p>
Formulation	Lyophilized from a 0.2 µm filtered solution of 10mM Tris, pH 8.0.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only