

# Recombinant Human Carbonic Anhydrase 14,CA14 (N-6His)

Catalog No: #AP75838

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP75838-1 10ug #AP75838-2 50ug #AP75838-3 500ug #AP75838-4 1mg

## Description

Product Name	Recombinant Human Carbonic Anhydrase 14,CA14 (N-6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9ULX7
Uniprot	Q9ULX7
GeneID	23632;
Calculated MW	32.8kDa
Target Sequence	MNHKVVHHHHHMGQHWTYEGPHGQDHPASYPECGNNAQSPIDIQTDSVTFDPDLPALQPHGYDQPGTEP LDLHNNGHTVQLSLPSTLYLGGLPRKYVAAQLHLHWGQKGSPPGGSEHQINSEATFAELHIVHYDSDSYDSLSE AAERPQGLAVLGILIEVGETKNIAYEHILSHLHEVRHKDQKTSVPPFNLRELLPKQLGQYFRYNGSLTTPPCYQS VLWTVFYRRSQISMEQLEKLGTLFSTEEEPSKLLVQNYRALQPLNQRMVFASFIQAGSSYTTGEMVDLQSR
Formulation	Supplied as a 0.2 µm filtered solution of 20mM Tris, 150mM NaCl, 10% Glycerol, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only