

## Recombinant Human Protein S100-A13,S100A13

Catalog No: #AP75841



Package Size: #AP75841-1 10ug #AP75841-2 50ug #AP75841-3 500ug #AP75841-4 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human Protein S100-A13,S100A13
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q99584
Uniprot	Q99584
GeneID	6284;
Calculated MW	11.3kDa
Target Sequence	AAEPLTELEESIETVVTFFTFARQEGRKDSLVSNEFKELVTQQLPHLLKDVGSLDEKMKSLDVNQDSELKFNE YWRLLIGELAKEIRKKKDLKIRKK
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB,150mM NaCl,pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only