

Recombinant Human Leukocyte Surface Antigen CD47,IAP,OA3 (C-Fc)



Catalog No: #AP75842

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75842-1 10ug #AP75842-2 50ug #AP75842-3 500ug #AP75842-4 1mg

Description

Product Name	Recombinant Human Leukocyte Surface Antigen CD47,IAP,OA3 (C-Fc)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q08722
Uniprot	Q08722
GeneID	961;
Calculated MW	40.8kDa
Target Sequence	QLLFNKTKSVEFTFCNDTVIPCFVTNMEAQNTTEVYVKWKFGRDIYTFDGalNKSTVPTDFSSAKIEVSQLL KGDASLKMdkSDAVSHTGNYTCEVTELTREGETIIElKYRVVSWFSPVDDIEGRMDEPKSCDKTHTCPPCPAP ELLGGPSVFLFPPKPKDTLMISRTPEVTCVVVDVSHEDPEVKFNWYVDGVEVHNAKTKPREEQYNSTYRVVS VLTVLHQDWLNGKEYKCKVSNKALPAPIEKTKAKGQPREPQVYTLPPSREEMTKNQVSLTCLVKGFYPSDIA VEWESNGQPENNYKTPPVLDSDGSFFLYSKLTVDKSRWQQGNVFCSCVMHEALHNYHTQKSLSLSPGK
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB150mM NaCl,pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only