

Recombinant Human Caspase-14,CASP14 (C-6His)

Catalog No: #AP75843



Package Size: #AP75843-1 10ug #AP75843-2 50ug #AP75843-3 500ug #AP75843-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Caspase-14,CASP14 (C-6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P31944
Uniprot	P31944
GeneID	23581;
Calculated MW	28.7kDa
Target Sequence	MSNPRSLEEEKYDMSGARLALILCVTKAREGSEEDLDALEHMFRQLRFESTMKRDPTAEQFQEELEKFQQAID SREDPVSCAFVVLMAHGREGFLKGEDGEMVKLENLFEALNNKNCQALRAKPKVYIIQACRGEQRDPGETVGG DEIVMVIKDSPQTIPTYTDALHVYSTVEGYIAYRHDQKQSCFIQTLVDVFTKRKGHILELLTEVTRRMAEAEVLQ EGKARKTNPEIQSTLRKRLYLQLEHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB150mM NaCl,pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only