

# Recombinant Human Ubiquitin-Conjugating Enzyme E2 R2,UBE2R2,UBC3B (N-6His)



Catalog No: #AP75846

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP75846-1 10ug #AP75846-2 50ug #AP75846-3 500ug #AP75846-4 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human Ubiquitin-Conjugating Enzyme E2 R2,UBE2R2,UBC3B (N-6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q712K3
Uniprot	Q712K3
GeneID	54926;
Calculated MW	29.3kDa
Target Sequence	MGSSHHHHHHSSGLVPRGSHMAQQMTSSQKALMLELKSLQEELVEGFRITLVDES DLYNWEVAIFGPPNTL YEGGYFKAHIKFPIDYPSPPTFRFLTKMWHPNYENG DVCISILHPPVDDPQSGELPSE RWNPTQNVRTILLS VISLLNEPNTFSPANVDASVMFRKWRDSKGDKEYAEIIRKQVSATKAEAEKDGKVPPTTLAEYCIKTKVPSND NSSDLLYDDLYDDDIDDEDEEEEDADCYDDDDSGNEES
Formulation	Supplied as a 0.2 µm filtered solution of 50mM HEPES, 150mM NaCl, pH 7.5.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only