

Recombinant Human Peroxiredoxin-3, PRDX3

Catalog No: #AP75848



Package Size: #AP75848-1 10ug #AP75848-2 50ug #AP75848-3 500ug #AP75848-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Peroxiredoxin-3, PRDX3
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P30048
Uniprot	P30048
GeneID	10935;
Calculated MW	21.6kDa
Target Sequence	MPAVTQHAPYFKGTAVVNGEFKDLSDDFKGGKYLVLFFYPLDFTFVCPTEIVAFSDKANEFHDVNCEVVAVSV DSHFSLAWINTPRKNGGLGHMNIALLSDLTKQISRDYGVLLLEGSLALRGLFIIDPNGVIKHLVNDLPVGRSV EETLRLVKAFQYVETHGEVCPANWTPDSPTIKPSPAASKEYFQKVNQ
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, pH8.0.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only