

# Recombinant Mouse C-C Motif Chemokine 21a,CCL21a,,6Ckine

Catalog No: #AP75849

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP75849-1 10ug #AP75849-2 50ug #AP75849-3 500ug #AP75849-4 1mg

## Description

Product Name	Recombinant Mouse C-C Motif Chemokine 21a,CCL21a,,6Ckine
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	P84444
Uniprot	P84444
GeneID	100504362;18829;
Calculated MW	12kDa
Target Sequence	SDGGGQDCCLKYSQKKIPYSIVRGYRKQEPSLGCPAILFSPRKHLSKPELCANPEEGWVQNLMRRLDQPPA PGKQSPGCRKNRGTSKSGKKGKGSKGCKRTEQTQPSRG
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB150mM NaCl,pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only