

Recombinant Human Ubiquitin, ISG15-Conjugating Enzyme E2 L6, UBE2L6, Ubch8 (C-6His)



Catalog No: #AP75851

Orders: order@signalwayantibody.com

Package Size: #AP75851-1 10ug #AP75851-2 50ug #AP75851-3 500ug #AP75851-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Ubiquitin, ISG15-Conjugating Enzyme E2 L6, UBE2L6, Ubch8 (C-6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	O14933
Uniprot	O14933
GeneID	9246;
Calculated MW	18.8kDa
Target Sequence	MMASMRVVKELEDLQKKPPPYLRNLSSDDANVLVWHALLLPDQPPYHLKAFNLRISFPPEYPFKPPMIKFTTKI YHPNVDENGQICLPIISSENWKPCTKTCQVLEALNVLVNRPNIREPLRMDLADLLTQNPFLFRKNAEEFTLRFG VDRPSLEHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 50mM HEPES, 100mM NaCl, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only