

Recombinant Human Tryptophan 5-Hydroxylase 2,TPH2 (N-6His)

Catalog No: #AP75852

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75852-1 10ug #AP75852-2 50ug #AP75852-3 500ug #AP75852-4 1mg

Description

Product Name	Recombinant Human Tryptophan 5-Hydroxylase 2,TPH2 (N-6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q14CC5
Uniprot	Q14CC5
Calculated MW	58.8kDa
Target Sequence	MNHKVVHHHHMMPAMMMFSSKYWARRGFLDSAVPEEHQLLGSSTASTFLKLNKPNKSGKNDKGNKGS KREAATESGKTAVVFLKNEVGGLVKALRLFQEKRVNMVHIESRKSRRRSSEVEIFVDCECGKTEFNELIQLLK FQTTIVTLNPPENIWTEEELEDPWFPRKISELDKCSHRVLMYGSELDADHPGFKDNVYRQRKYFVDVAMG YKYGQPIPRVEYTEEETKTWGVVFRELSKLYPTHACREYLKNFPLLTKYCGYREDNVPQLEDVSMFLKERSGF TVRPVAGYLSRDFLAGLAYRVFHCTQYIRHGSDPLYTPEPDTHELLGHVPLLADPKFAQFSQEIGLASLGAS DEDVQKLATCYFFTIEFGLCKQEGQLRAYGAGLLSSIGELKHALSDKACVKAFDPKTTCLQECLITTFQEAYFV SESFEEAKEKMRDFAKSITRPFVYFNPYQTSIEILKDTRSIENVVQDLRSDLNTVCDALNKMNQYLGIVDLQSR
Formulation	Supplied as a 0.2 µm filtered solution of 20mM PB150mM NaCl,pH7.4.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only