

Recombinant Human Ubiquitin-Conjugating Enzyme E2 Z,UBE2Z (N-6His)



Catalog No: #AP75854

Orders: order@signalwayantibody.com

Package Size: #AP75854-1 10ug #AP75854-2 50ug #AP75854-3 500ug #AP75854-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Ubiquitin-Conjugating Enzyme E2 Z,UBE2Z (N-6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9H832-2
Uniprot	Q9H832
GeneID	65264;
Calculated MW	30.2kDa
Target Sequence	MGSSHHHHHSSGLVPRGSHMSIYKEPPPGMFVVPDVTDMTKIHALITGPFDTPTYEGGFLLVFRCPDPIH PPRVKLMTTGNNTVRFNPNFYRNGKVCLSILGTWTGPAWSPAQSISVLSIQSLMTENPYHNEPGFEQERHP GDSKNYNECIRHETIRVAVCDMMEGKPCPEPLRGVMEKSFLEYDYDFYEVACKDRLHLQGQTMQDPFGEKR GHFDYQSLLMRLGLIRQKVLRLHNENAEMDSDSSSSGTETDLHGSLRV
Formulation	Supplied as a 0.2 µm filtered solution of 50mM HEPES, 50mM NaCl, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only