

Axl (Phospho-Tyr698) Polyclonal Antibody

Catalog No: #12335



Package Size: #12335-1 50ul #12335-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Axl (Phospho-Tyr698) Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC
Species Reactivity	Hu Ms Rt
Specificity	Phospho-Axl (Y698) Polyclonal Antibody detects endogenous levels of Axl protein only when phosphorylated at Y698.
Immunogen Description	Synthesized peptide derived from human Axl around the phosphorylation site of Y698.
Target Name	Axl
Modification	Phospho
Other Names	AXL; UFO; Tyrosine-protein kinase receptor UFO; AXL oncogene
Accession No.	Swiss-Prot: P30530NCBI Gene ID: 558
Uniprot	P30530
GeneID	558;
Target Species	human
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

WB 1:500-1:2000

IHC 1:100-1:300

Note: This product is for in vitro research use only