

Recombinant Human Signal Transducer and Activator of Transcription 5B, STAT5B (C-6His)



Catalog No: #AP75857

Orders: order@signalwayantibody.com

Package Size: #AP75857-1 10ug #AP75857-2 50ug #AP75857-3 500ug #AP75857-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Signal Transducer and Activator of Transcription 5B, STAT5B (C-6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P51692
Uniprot	P51692
GeneID	6777;
Calculated MW	38.4kDa
Target Sequence	MAVWIQAQQLQGEALHQMQUALYQGHFPIEVRHYLSQWIESQAWDSVDLDNPPQENIKATQLLEGLVQELQKKA EHQVGEDGFLLKIKLGHYATQLQNTYDRCPMELVRCIRHILYNEQRLVREANNNGSSPAGSLADAMSQKHLQIN QTFEELRLVTQDTENELKKLQQTQEYFIIQYQESLRQAQFGPLAQLSPQERLSRETALQQKQVSLEAWLQRE AQTLQQYRVELAEKHQKTLQLLRKQQTIIIDDELIQWKRRQQLAGNGGPPEGSLDVLQSWCEKLAEIWQNRQ QIRRAEHLCCQLPIPGPVEEMLAEVNATITLHHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of PBS, 50% Glycerol, 1mM DTT, pH 7.4.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only