

Recombinant Human Ubiquitin-Like Protein ISG15,ISG15 (C-6His)

Catalog No: #AP75859

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75859-1 10ug #AP75859-2 50ug #AP75859-3 500ug #AP75859-4 1mg

Description

Product Name	Recombinant Human Ubiquitin-Like Protein ISG15,ISG15 (C-6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P05161
Uniprot	P05161
GeneID	9636;
Calculated MW	18.2kDa
Target Sequence	<p> GWDLTVKMLAGNEFQVSLSSSMSVSELKAQITQKIGVHAFQQLAVHPSGVALQDRVPLASQGLGPGSTVLL VVDKCDEPLNILVRNNKGRSSTYEVRLTQTVAHLKQQVSGLEGVQDDLFWLTFEGKPLEDQLPLGEYGLKPLS TVFMNLRRLRGGLEHHHHHHH </p>
Formulation	Supplied as a 0.2 µm filtered solution of 50mM HEPES, 100mM NaCl, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only