

Recombinant Human α -Taxilin, TXLNA (N, C-6His)

Catalog No: #AP75860



Package Size: #AP75860-1 10ug #AP75860-2 50ug #AP75860-3 500ug #AP75860-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human α -Taxilin, TXLNA (N, C-6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P40222
Uniprot	P40222
GeneID	200081;
Calculated MW	20.4kDa
Target Sequence	MGSSHHHHHHSSGLVPRGSHMKNQDKKNGAAKQSNPKSSPGQPEAGPEGAQERPSQAAPAVEAEGPGSS QAPRKPEGAQARTAQSGALRDVSEELSRQLEDILSTYCVDNQGGPGEDGAQGEPAEPEDA EKSRTYVARN GEPEPTPVNNGEKEPSKGDPNTEEIRQSDEVGDRDHRPQEKLEHHHHHH
Formulation	Lyophilized from a 0.2 μ m filtered solution of 20mM PB, 150mM NaCl, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only