

Recombinant Human Fibroblast Growth Factor 17,FGF-17

Catalog No: #AP75861

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75861-1 10ug #AP75861-2 50ug #AP75861-3 500ug #AP75861-4 1mg

Description

Product Name	Recombinant Human Fibroblast Growth Factor 17,FGF-17
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	O60258
Uniprot	O60258
GeneID	8822;
Calculated MW	22.6kDa
Target Sequence	MTQGENHPSNFNQYVRDQGAMTDQLSRRQIREYQLYSRTSGKHVQVTGRRISATAEDGNKFAKLIVETDTF GSRVRIKGAESEKYICMNRKGLIGKPSGKSKDCVFTEIVLENNYAFQNRHEGWMAFTRQGRPRQASRS RQNQREAHFIKRLYQGQLPFPNHAEKQKQFEFVGSAPTRRTKRTRRPQPLT
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB,150mM NaCl,pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only