

Recombinant Human Parvulin-14,PIN4 (N-6His)

Catalog No: #AP75862



Package Size: #AP75862-1 10ug #AP75862-2 50ug #AP75862-3 500ug #AP75862-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Parvulin-14,PIN4 (N-6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9Y237-2
Uniprot	Q9Y237
GeneID	5303;
Calculated MW	18.8kDa
Target Sequence	MGSSHHHHHHSSGLVPRGSHMPMAGLLKGLVRQLEQFRVQQASKMPPKKGKSGSGKAGKGAASGSDSA DKKAQGPKGGGNAVKVRHILCEKHGKIMEAMEKLKSGMRFNVEAAQYSEDKARQGGDLGWMTRGSMVGPFF QEAAFALPVSGMDKPVFTDPPVKTKFGYHIIMVEGRK
Formulation	Supplied as a 0.2 µm filtered solution of 20mM PB,150mM NaCl,pH7.4.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only