

# Recombinant Human Ubiquitin-Conjugating Enzyme E2 K,UBE2K (N-6His, SUMO tag)



Catalog No: #AP75863

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP75863-1 10ug #AP75863-2 50ug #AP75863-3 500ug #AP75863-4 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human Ubiquitin-Conjugating Enzyme E2 K,UBE2K (N-6His, SUMO tag)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P61086
Uniprot	P61086
GeneID	3093;
Calculated MW	34.5kDa
Target Sequence	MGHHHHHHGSDSEVNQEAKPEVKPEVKPETHINLKVSDGSSEIFFKIKKTTPLRRLMEAFKRQGGKEMDSLRF LYDGIRIQADQAPEDLDMEDNDIIEAHREQIGGMANIAVQRIKREFKEVLKSEETSKNQIKVDLVDFTELGEI AGPPDTPYEGGRYQLEIKIPETYPFNPPKVRFITKIWHPNISSVTGAICLDILKDQWAAAMTLRTVLLSLQALLAA AEPDDPQDAVVANQYKQNPPEMFKQTARLWAHVYAGAPVSSPEYTKKIENLCAMGFDRNAVIVALSSKSWDV ETATELLLSN
Formulation	Supplied as a 0.2 µm filtered solution of 20mM PB,150mM NaCl,pH7.4.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only