

Recombinant Human Syntenin-1,SDCBP,SYCL,MDA9 (C-6His)



Catalog No: #AP75864

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75864-1 10ug #AP75864-2 50ug #AP75864-3 500ug #AP75864-4 1mg

Description

| | |
|-----------------|---|
| Product Name | Recombinant Human Syntenin-1,SDCBP,SYCL,MDA9 (C-6His) |
| Host Species | E. coli |
| Purification | Greater than 95% as determined by reducing SDS-PAGE. |
| Specificity | Human |
| Accession No. | O00560 |
| Uniprot | O00560 |
| GeneID | 6386; |
| Calculated MW | 33.5kDa |
| Target Sequence | SLYPSLEDLKVVDKVIQAQTAFSANPANPAILSEASAPIPHDGNLYPRLYPELSQYMGLSLNEEEIRANVAVVSGA PLQGQLVARPSSINYMVAPVTGNDVGIIRRAEIKQGIREVILCKDQDQDGKIGLRLKSIDNGIFVQLVQANSPASLVG LRFQDQVLQINGENCAGWSSDKAHKVLKQAFGEKITMTIRDRPFERTITMHKDSTGHVGFIFKNGKITSIVKDS SAARNGLLTEHNICEINGQNVIGLKDSQIADILSTSGTVVTITIMPAFIFEHIIKRMAPSIMKSLMDHTIPEVLEHHH HHH |
| Formulation | Supplied as a 0.2 µm filtered solution of 20mM PB,150mM NaCl,pH7.4. |
| Storage | Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles. |

Note: This product is for in vitro research use only