

Recombinant Human NEDD8-Conjugating Enzyme UBC12,UBE2M



Catalog No: #AP75867

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75867-1 10ug #AP75867-2 50ug #AP75867-3 500ug #AP75867-4 1mg

Description

Product Name	Recombinant Human NEDD8-Conjugating Enzyme UBC12,UBE2M
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P61081
Uniprot	P61081
GeneID	9040;
Calculated MW	20.9kDa
Target Sequence	MIKLFSLKQQKKEEESAGGKGGSSKKASAAQLRIQKDINELNLPKTCDISFSDPDDLNFKLVICPDEGFYKSGK FVFSFKVGGYPHDPPKVKCETMVYHPNIDLEGNVCLNILREDWKPVLINSIIYGLQYLFLEPNPEDPLNKEAA EVLQNNRRRLFQNVQRSMRGGYIGSTYFERCLK
Formulation	Supplied as a 0.2 µm filtered solution of 50mM HEPES,pH 7.5,2mM DTT,150mM NaCl,10%glycerol .
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only