

# Recombinant Mouse Tumor Necrosis Factor Receptor I, TNFRSF1A, CD120a



Catalog No: #AP75868

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP75868-1 10ug #AP75868-2 50ug #AP75868-3 500ug #AP75868-4 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Mouse Tumor Necrosis Factor Receptor I, TNFRSF1A, CD120a
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	P25118
Uniprot	P25118
GeneID	21937;
Calculated MW	21.2kDa
Target Sequence	MIHPSGVTGLVPSLGDREKRDSLCPQGKYVHSKNNSICCTKCHKGTYLVSDCPSPGRDTVCRECEKGTFTAS QNYLRQCLSKCTCRKEMSQVEISPCQADKDTVCGCKENQFQRYLSETHFQCVDCCSPCFNGTVTIPCKETQNT VCNCHAGFFLRESECVPCSHCKKNEECMKLCLPPPLANVTNPQDSGTA
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB150mM NaCl, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only