

## Recombinant Human Resistin,RETN,ADSF (C-6His)

Catalog No: #AP75869



Package Size: #AP75869-1 10ug #AP75869-2 50ug #AP75869-3 500ug #AP75869-4 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human Resistin,RETN,ADSF (C-6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9HD89
Uniprot	Q9HD89
GeneID	56729;
Calculated MW	10.6kDa
Target Sequence	KTLCSMEEAINERIQEVAGSLIFRAISSIGLEQSVTSRGLATCPRGFAVTGCTCGSACGSWDVRAETTCHC QCAGMDWTGARCCRVQPLEHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM Acetic acid,pH 3.0.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only