

Recombinant Human Ankyrin Repeat and SOCS Box Protein 13,ASB13 (N-6His)



Catalog No: #AP75870

Orders: order@signalwayantibody.com

Package Size: #AP75870-1 10ug #AP75870-2 50ug #AP75870-3 500ug #AP75870-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Ankyrin Repeat and SOCS Box Protein 13,ASB13 (N-6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q8WXK3
Uniprot	Q8WXK3
GeneID	79754;
Calculated MW	31.4kDa
Target Sequence	MNHKVVHHHHHMEPRAADGCFLGDVGFVVVERTPVHEAAQRGESLQLQQLIESGACVNQVTVDSTPLHAAS LQQQARCVQLLLAAGAQVDARNIDGSTPLCDACASGSIECVKLLLSYGAKVNPPLYTASPLHEACMSGSSSECV RLLIDVGANLEAHDCHFGLPLHVACAREHLDCVKVLLNAGANVNAAKLHETALHHAAKVKNVDLIEMLIIEFGGN IYARDNRGKKPSDYTWSSAPAKCFEYYEKTPLTLSQLCRVNLKATGVRGLEKIAKLNIPRLIDYLSYN
Formulation	Supplied as a 0.2 µm filtered solution of PBS, pH 7.4, 4M Urea.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only