

Recombinant Human Nogo-A,Reticulon-4,RTN4 (N-GST)

Catalog No: #AP75873

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75873-1 10ug #AP75873-2 50ug #AP75873-3 500ug #AP75873-4 1mg

Description

Product Name	Recombinant Human Nogo-A,Reticulon-4,RTN4 (N-GST)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9NQC3-3
Uniprot	Q9NQC3
GeneID	57142;
Calculated MW	45.5kDa
Target Sequence	MSPILGYWKIKGLVQPTRLLEYLEEKYEEHLYERDEGDKWRNKKFELGLEFPNLPYYIDGDVKLTQSMAIIRYI ADKHNMLGGCPKERAEISMLEGAVLDIRYGVSR IAYS KDFETLKVDFLSKLP EMLKMFEDRLCHKTYLNGDHV THPDFMLYDALDVVLYMDPMCLDAFPKLVCFKKRIEAI PQIDKYLKSSKYIAWPLQGWQATFGGGDHPPKSDL VPRGSMEDLDQSPLVSSSDSPRPQPAFKYQFVREPEDEEEEEEEEEDEDEDLEELEVLERKPAAGLSAAP VPTAPAAGAPLMDFGNDVFPAPRGPLPAAPPVAPERQPSWDPSPVSSTVPAPSPLSAAAVSPSKLPEDDEP PARPPPPPPASVSPQAEPVWTPPAPAPAAPPSTPAAPKRRGSSGSV
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB150mM NaCl,pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only