

# Recombinant Human Chloride Intracellular Channel Protein 3, CLIC3 (C-6His)

Catalog No: #AP75875

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP75875-1 10ug #AP75875-2 50ug #AP75875-3 500ug #AP75875-4 1mg

## Description

Product Name	Recombinant Human Chloride Intracellular Channel Protein 3, CLIC3 (C-6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	O95833
Uniprot	O95833
GeneID	9022;
Calculated MW	27.7kDa
Target Sequence	MAETKLQLFVKASEDGESVGHCPSCQRLFMVLLKGVPFLLTTVDTRRSPDVLKDFAPGSQLPILLYDSDAKT DTLQIEDFLEETLGPPDFPSLAPRYRESNTAGNDVFHKFSAFIKNPVPAQDEALYQQLLRALARLDSYLRAPLE HELAGEPQLRESRRRFLDGDRLLADCSLLPKLHIVDTVCAHFRQAPIPAELRGVRRYLD SAMQEKEFKYTCP HSAEILAAAYRPAVHPRLEHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 10mM Tris, 0.1% Triton100, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only