

Recombinant Human Endoglin,CD105 (N-Trx, 6His)

Catalog No: #AP75878



Package Size: #AP75878-1 10ug #AP75878-2 50ug #AP75878-3 500ug #AP75878-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Endoglin,CD105 (N-Trx, 6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P17813
Uniprot	P17813
GeneID	2022;
Calculated MW	33.6kDa
Target Sequence	MSDKIIHLTDDSFDTDLKADGAILVDFWAEWCGPCKMIAPILDEIADEYQGKLTVAKLNIDQNPGTAPKYGIRGI PTLLLFKNGEVAATKVGALSKGQLKEFLDANLAGSGSGHMHSHHHSSGLVPRGSGMKETAATAAKFERQHMD SPDLGTDDDDKAMETVHCDLQPVGPERDEVITYTTSQVSKGCVAQAPNAILEVHVLFLFEPFGPSQLELTLQAS KQNGTWPREVLLVLSVNSSVFLHLQALGIPLHLAYNSSLVTFQEPGPNTELPSPFKTQILEWAAERGPITSA AELNDPQSILLRLGQAQ
Formulation	Supplied as a 0.2 µm filtered solution of 20mM PB150mM NaCl,pH7.4.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only