

Recombinant Human Calreticulin-3,CALR3,CRT2

Catalog No: #AP75879



Package Size: #AP75879-1 10ug #AP75879-2 50ug #AP75879-3 500ug #AP75879-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Calreticulin-3,CALR3,CRT2
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q96L12
Uniprot	Q96L12
GeneID	125972;
Calculated MW	42.9kDa
Target Sequence	MTVYFQEEFLDGEHWRNRWLQSTNDSRFGHFRLSSGKFGYGHKEKDKGLQTTQNGRFYAISARFKPFSNKGK TLVIQYTVKHEQKMDCGGGYIKVFPADIDQKNLNGKSQYYIMFGPDICGFDIKKVHVILHFKNKYHENKKLIRCK VDGFTHLYTLILRPDLSYDVKIDGQSIESGSIEYDWNLTSLKKETSPAESKDWEQTKDNKAQDWEKHFLDASTS KQSDWNGDLGDWPAPMLQPPYQDGLKPEGIHKDVLHRKMKNTDYLTQYDLSEFENIGAIGLELWQVRS GTIFDNFLITDDEEYADNFGKATWGETKGPEREMDAIQAKEEMKKAREEEEEELLSGKINRHEHYFNQFHRRN EL
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB150mM NaCl,pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only