

Recombinant Human Guanidinoacetate N-Methyltransferase, GAMT (N, C-6His)



Catalog No: #AP75880

Orders: order@signalwayantibody.com

Package Size: #AP75880-1 10ug #AP75880-2 50ug #AP75880-3 500ug #AP75880-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Guanidinoacetate N-Methyltransferase, GAMT (N, C-6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q14353
Uniprot	Q14353
GeneID	2593;
Calculated MW	29.5kDa
Target Sequence	MGSSHHHHHHSSGLVPRGSHMSAPSATPIFAPGENCSPAWGGAAPAAAYDAADTHLRILGKPVMERWETPYMH ALAAAASSKGGRVLEVGFMAIAASKVQEAPIDEHWIIECNDGVFQRLRDWAPRQTHKVIPLKGLWEDVPTL PDGHFDGILYDITYPLSEETWHTHQFNFIKNHAFRLKPGGVLTYCNLTSWGELMKSKEYSDITIMFEETQVPALL EAGFRRENIRTEVMALVPPADCRYYAFPQMITPLVTKGLEHHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 1mMDTT, pH 8.0 .
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only