

Recombinant Human Sterol O-Acyltransferase 2,SOAT2,ACAT2 (N-Trx, 6His)



Catalog No: #AP75894

Orders: order@signalwayantibody.com

Package Size: #AP75894-1 10ug #AP75894-2 50ug #AP75894-3 500ug #AP75894-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Sterol O-Acyltransferase 2,SOAT2,ACAT2 (N-Trx, 6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	O75908
Uniprot	O75908
GeneID	8435;
Calculated MW	59.1kDa
Target Sequence	MSDKIIHLTDDSFDTDLKADGAILVDFWAEWCGPCKMIAPILDEIADEYQGKLTVAKLNIDQNPGTAPKYGIRGI PTLLLFKNGEVAATKVGALSKGQLKEFLDANLAGSGSGHMHSHHHSSGLVPRGSGMKETAALKFERQHMD SPDLGTDDDDKAMADIGSMNAGSDPVVIVSAARTIIGSFNGALAAVPVQDLGSTVIKEVLKRATVAPEDVSEVIF GHVLAAGCGQNPVRQASVAGIPYSVPAWSCQMIGSGLKAVCLAVQSIGIGDSSIVVAGGMENMSKAPHLA YLRTGVKIGEMPLTDSILCDGLTDAFHNMGITAENVAKKWQVSREDQDKVAVLSQNRTEAQAQKAGHFDKEI VPVLVSTRRLIEVKTDEFPRHGSNIEAMSKLPYFLDTGTGTVPANASGINDGAAAVVLMKKSEADKRGLTP LARIVSWSQVGVPEPSIMGIGPIPAIKQAVTKAGWSLEDVDIFEINEAFAAVSAAIVKELGLNPEKVNIEGGAIALG HPLGASGCRILVTLHHTLERMGRSRGVAALCIGGGMGIAMCVQRE
Formulation	Supplied as a 0.2 µm filtered solution of 20mM PB150mM NaCl,pH7.4.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only