

Recombinant Human Nucleolar Protein Family A Member 2,NHP2,NOLA2 (N-6His)

Catalog No: #AP75896

Orders: order@signalwayantibody.com

Package Size: #AP75896-1 10ug #AP75896-2 50ug #AP75896-3 500ug #AP75896-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Nucleolar Protein Family A Member 2,NHP2,NOLA2 (N-6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9NX24
Uniprot	Q9NX24
GeneID	55651;
Calculated MW	19.3kDa
Target Sequence	MGSSHHHHHHSSGLVPRGSHMTKIKADPDGPEAQAEACSGERTYQELLVNQNPIAQPLASRRLTRKLYKCIK KAVKQKQIRRGVKEVQKFVNKGEKIMVLAGDTLPIEVYCHLPVMCEDRNLPLYVYIPSKTDLGAAAGSKRPTC VIMVKPHEEYQEAYDECLEEVQSLPLPL
Formulation	Supplied as a 0.2 µm filtered solution of 20mM Tris,100mM NaCl,1mM DTT, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only