

Recombinant Human Peptidase D,PEPD

Catalog No: #AP75899



Package Size: #AP75899-1 10ug #AP75899-2 50ug #AP75899-3 500ug #AP75899-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Peptidase D,PEPD
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P12955
Uniprot	P12955
GeneID	5184;
Calculated MW	54.5kDa
Target Sequence	AAATGPSFWLGNETLKVPLALFALNRQRLCERLRKNPAVQAGSIVVLQGGREETQRYCTDTGVLFRRQESFFHW AFGVTEPGCYGVIDVDTGKSTLFPRLPASHATWWMGKIHSKEHFKEKYAVDDVQYVDEIASVLTSQKPSVLLTL RGVNTDSGSVCREASFDGISKFEVNTILHPEIVECRVFKTDMELEVLRYTNKISSEAHREVMKAVKVGKEY ELESLEFYCYSRGGMRHSSYTCICGSGENSAVLHYGHAGAPNDRTIQNGDMCLFDMGGEYYCFASDITCSF PANGKFTADQKAVYEAVLRSSRAVMGAMKPGVWWPDMHRLADRIHLEELAHMGILSGSVDAMVQAHLGAVF MPHGLGHFLGIDVHDVGGYPEGVERIDEPGLRSLRTARHLQPGMVLVEPGIYFIDHLLDEALADPARASFFN REVLQRFRGFGGVRIEEDVVVTDSGIELLTCVPRTVEEIEACMAGCDKAFTPFSGPK
Formulation	Supplied as a 0.2 µm filtered solution of 20mM PB150mM NaCl,pH7.4.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only