

Recombinant Human Interleukin-1 β , IL-1 β

Catalog No: #AP75904



Package Size: #AP75904-1 10ug #AP75904-2 50ug #AP75904-3 500ug #AP75904-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Interleukin-1 β , IL-1 β
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P01584
Uniprot	P01584
GeneID	3553;
Calculated MW	17.5kDa
Target Sequence	MAPVRSLNCTLRDSQQKSLVMSGPYELKALHLQGQDMEQQVVFMSFVQGEESNDKIPVALGLKEKNLYLSC VLKDDKPTLQLESVDPKNYPKKKMEKRFVFNKIEINNKLFEFSAQFPNWIYSTQAENMPVFLGGTKGGQDIT DFTMQFVSS
Formulation	Lyophilized from a 0.2 μ m filtered solution of 20mM PB, 150mM NaCl, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only