

Recombinant Human Eukaryotic Translation Initiation Factor 5A-2, EIF5A2 (N-6His)



Catalog No: #AP75906

Orders: order@signalwayantibody.com

Package Size: #AP75906-1 10ug #AP75906-2 50ug #AP75906-3 500ug #AP75906-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Eukaryotic Translation Initiation Factor 5A-2, EIF5A2 (N-6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9GZV4
Uniprot	Q9GZV4
GeneID	56648;
Calculated MW	18.9kDa
Target Sequence	MGSSHHHHHHSSGLVPRGSHMADEIDFTTGDAGASSTYPMQCSALRKNQFVVLKGRPCKIVEMSTSKTGKH GHAKVHLVGIDIFTGKKYEDICPSTHNMDVPMNIKRNDYQLICIQDGYLSLLTETGEVREDLKLPEGELGKEIEGK YNAGEDVQVSVMCAMSEEYAVAIKPKK
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM Tris, 150mM NaCl, pH 8.0.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only