

# Recombinant Cavia porcellus Interleukin-1 Beta, IL-1 beta (C-6His)

Catalog No: #AP75907

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP75907-1 10ug #AP75907-2 50ug #AP75907-3 500ug #AP75907-4 1mg

## Description

Product Name	Recombinant Cavia porcellus Interleukin-1 Beta, IL-1 beta (C-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Cavia porcellus
Accession No.	Q9WVG1
Uniprot	Q9WVG1
GeneID	100135556;
Calculated MW	18.2kDa
Target Sequence	MTPVPSRNCTLHDIQHKRLVLSDPCELKALHLNGDNLNRQVVFMSFVQGERSDNKMPVALGLKGKNLYLSC VMKDGKPVQLQLESVDGKQYPKKKMEKRFVFNKITSKSTVEFESAQFPNWIYSTSQAEHKPVFLGNNGQDIID FKLELVSSHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only