

# Recombinant Human Microtubule-associated Proteins 1A,1B Light Chain 3B,MAP1LC3B



Catalog No: #AP75908

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP75908-1 10ug #AP75908-2 50ug #AP75908-3 500ug #AP75908-4 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human Microtubule-associated Proteins 1A,1B Light Chain 3B,MAP1LC3B
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9GZQ8
Uniprot	Q9GZQ8
GeneID	81631;
Calculated MW	14.8kDa
Target Sequence	GSMPSEKTFKQRRTFEQRVEDVRLIREQHPTKIPVIERKGEKQLPVLDKTKFLVPDHVNMSELIKIIRRLQLN ANQAFFLLVNGHSMVSVSTPISEVYESEKDEDGFLYMVYASQETFGMKLSV
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM Tris, 150mM NaCl, 2mM DTT, pH8.0.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only